## Product Anatomy and Associated Parts FLT-KB##



Stealth Barbed Watering System





FLT-030 DIN Vent Adapter for Stealth Valves, (24 pack)

FLT-068 Twist-Lock Male Input Coupling for Stealth Watering Systems (6 pack)



FLT-093 Tubing Cutter for Speed Tubing (Water Injector Spider & Stealth Watering Systems)



FLT-099 Standard Input Coupling and Filter Assembly for Stealth Barbed Watering System, 10 pack



FLT-108 Tubing, 100', Black PVC Tubing, Barbed Stealth



FLT-104 Input Filter used in Stealth Barbed Input Assembly, with collars (6 pack)



FLT-105 End Caps for Stealth Barbed Watering Systems, 50 pack



FLT-BLD1 Stealth Barbed Valve, Single, with DIN Adapter

FLT-BL1 Barbed Stealth Float, Loose, Black



INJ-198 Tubing, 25 feet 1/2" OD






INJ-601 Collar Installation Tool



INJ-608 Collars for Water Injectors & Stealth Systems (25 Pack)



FLT-117 Tubing, 200', Black PVC, Stealth Barbed



Philadelphia Scientific<sup>®</sup> industrial battery innovation

Philadelphia Scientific LLC 207 Progress Drive Montgomeryville, PA 18936 USA P: 215.616.0390 F: 215.616.0500 www.phlsci.com

DOC0430 4 / 2022

© 2022 Philadelphia Scientific LLC All Rights Reserved. Philadelphia Scientific and the PS logo are registered trademarks of Philadelphia Scientific LLC. Stealth is a trademark of Philadelphia Scientific LLC. E&O.E.